# VACUUM MATTRESS 4 872.07

CE

Regolamento (UE)  
2017/745

IT

[WWW.KONG.IT](http://WWW.KONG.IT)

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Declares under its sole responsibility that the device:

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Intended use:

Vacuum mattress for patient immobilization

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# INDEX

<b>1 - SYMBOLS AND SUPPORT</b>	<b>4</b>
1.1 Symbols	4
1.2 Support	4
<b>2 - GENERAL INFORMATION</b>	<b>5</b>
<b>3 - TECHNICAL FEATURES</b>	<b>6</b>
3.1 Parts terminology and materials	6
3.2 Dimension and capacity	8
3.3 Compatible devices, optional components and spare parts	8
<b>4 - SPECIFIC INFORMATION</b>	<b>10</b>
4.1 Applications	10
4.2 Preparation	10
4.3 Immobilising the patient	11
4.4 Trasporting the patient	12
4.5 Storage of the device	12
<b>5 - MAINTENANCE AND REPAIR</b>	<b>13</b>
5.1 General Information	13
5.2 Maintenance	13
5.3 Repair	13
<b>6 - STORAGE</b>	<b>14</b>
<b>7 - CHECKS, INSPECTIONS AND REVISIONS</b>	<b>15</b>
7.1 Checks	15
7.2 Inspections	15
7.3 Revisions	15
<b>8 - DURATION AND WARRANTY</b>	<b>16</b>
8.1 Duration	16
8.2 Disposal	16
8.3 Warranty	16
8.4 Law obligations	16
<b>9 - LABELLING AND SYMBOLS</b>	<b>17</b>
9.1 Labelling	17
9.2 Separable components label	17
9.3 Symbols	17
<b>10 - DOCUMENTS</b>	<b>18</b>
10.1 Inspections and servicing register	18
10.2 Maintenance and repair register	19
10.3 EU Declaration of conformity	20
<b>11 - LIST OF APPLIED AND REFERENCE STANDARDS</b>	<b>21</b>
11.1 Applied standards	21
11.2 Standards used as reference	22

# CHAPTER

# 1

## SYMBOLS AND SUPPORT

### 1.1 SYMBOLS

For the sake of comfort and clarity while reading this manual, below are the symbols used for handling important warnings for a proper and safe use of the device.

#### **REQUIREMENT FOR PROPER USE**

It identifies the presence of information for proper use of the device.

---

#### **INFORMATION REQUIREMENT**

It identifies the presence of useful and general information which guides the reader towards a conscious use of the device and/or performance of actions.

---

It identifies that the Medical Device is manufactured, designed and produced in compliance with the provisions of the General Safety and Performance Requirements of Regulation (EU) 2017/745 (Class I medical device, in compliance with classification rule 1 as indicated in the attachment VIII).

---

### 1.2 SUPPORT

For information please contact Kong Customer Support Service by:

- telephone +39 0341 630506
- fax +39 0341 641550
- email: [safetycare@kong.it](mailto:safetycare@kong.it)

or write to KONG S.p.A. – Via XXV Aprile, 4 – 23804 Monte Marenzo LC - ITALY.

To facilitate support operations, please always communicate or state the serial number (SN) indicated on the label fixed to the Medical Device.

# CHAPTER 2

## 2 GENERAL INFORMATION

The information provided by the manufacturer (hereinafter referred to as information) must be read and clearly understood by the user prior to using the Medical Device (MD hereafter). The information regards the description of the features, performance, assembly, disassembly, maintenance, preservation, disinfection, etc. of the device. Even though the information offers tips on use, this information shall not be considered as a user manual under actual conditions of use.

### WARNINGS AND LIMITATIONS OF USE:

- this device shall be strictly used by people who are physically fit, trained (instructed and taught) to use the device and with specific experience regarding moving the patient or, during training activities, by people under direct supervision of the trainers/ supervisors who guarantee the safety thereof,
- do not use the device before fully reading and understanding this user manual,
- strictly follow the manufacturer's information, improper use of the device is hazardous,
- modifying and/or repairing the device is strictly forbidden,
- all checks described in chapter 7 must be carried out prior to and after using the device. In case of any doubt on the efficiency of the device, the user must replace it immediately,
- non-compliant use, deformations, falls, wear, chemical contamination, exposing textile/plastic components/devices to temperature below -30°C or above

- +50°C and metal components/devices to temperatures exceeding 100°C are some examples of causes that can reduce, limit and end the life of the device,
- Prior to any rescue operation, be keen not to exceed the capacity indicated in paragraph 3.3,
- in order to reduce risks of exposure to / transmission of infectious diseases, clean and disinfect the device as indicated in chapter 5,
- improper use of the patient immobilisation systems can jeopardise the safety of the patient,
- always check the compatibility of the devices used alongside the device by consulting the relative manufacturer's information,
- use of spare parts or optional components different from the ones indicated in paragraph 3.4 can be hazardous,
- do not expose the device to sources of heat and at contact with chemical substances. Reduce exposure to direct sunlight as much as possible. At low temperatures and in the presence of humidity, formation of ice could reduce flexibility and increase cutting and abrasion-related risks on textile and synthetic devices,
- Report to the manufacturer and the competent authority of the member state where the user and / or patient is established any serious incident occurring in relation to the MD.

All our devices are tested/ checked piece by piece in compliance with the procedures laid down by the Quality System certified in accordance with the UNI EN ISO 9001 standard. Laboratory tests, testing, information and standards do not always reproduce the practical result. Thus, the results obtained under the actual conditions of use of the device in the natural environment may differ, even considerably at times. The best information lies in the continuous practical use under the supervision of skilled/ expert/qualified people.

**Warning: not suitable for use in an ATEX environment (Directive 94/9/CE)**

# CHAPTER 3

## TECHNICAL FEATURES

### 3.1 TERMINOLOGY AND MATERIALS OF THE PARTS

- A - PVC load bearing structure,
- B - Vacuum chambers containing expanded polyester spheres,
- C - Polypropylene locking straps,
- D - Polypropylene carrying handles,
- E - Polypropylene head immobilizer straps.



Fig.1

- F - PVC connection loops,
- F1 - PVC numbered connection loops,
- G - Ventral band in polypropylene,
- H - Openings for the passage of the straps of compatible devices,
- I - Valve,
- L - Aluminum alloy suction pump,
- M - Nylon bag.


Fig.1

## 3.2 DIMENSIONS AND CAPACITY

Dimensions:	185x70 cm
Dimensions (stored in the bag):	60x45x20 cm
Overall weight of the parts (fig.1):	6 kg
Load capacity:	150 kg
Carrying handles:	8 (4 on each side)
Independent vacuum chambers:	7, containing expanded polyester spheres
Radiotransparency:	≈ 920 U.H. Tested through radiographic and CT tests
Use/storage conditions:	Temperature: -20 / +60 C Relative humidity: 10 - 90 %

## 3.3 COMPATIBLE DEVICES, OPTIONAL COMPONENTS AND SPARE PARTS

Carefully read the pairing instructions in the manual for compatible devices and optional components.

### 3.3.1 Compatible devices

- 855.00 EVEREST BAG (using the VACUUM EVEREST KIT - 855.15)
- 870.04 LECCO 2.0 STRETCHER
- 880.01 911 CANYON STRETCHER
- 880.02 911 NET FULL STRETCHER
- 880.03 911 SHELL STRETCHER
- 890.01 ROLLY STRETCHER

### 3.3.2 Optional components

- 855.15 VACUUM EVEREST KIT

The kit allows you to connect the VACUUM MATTRESS 4 to the EVEREST BAG.


Fig.2

### 3.3.3 Spare parts

L - Suction pump,

M - Nylon bag.

# CHAPTER 4

## SPECIFIC INFORMATION

### 4.1 APPLICATION

The “VACUUM MATTRESS 4” is a medical device suitable for the immobilization of patients and intended for hospital use, for health emergencies, particularly in cases of long transport or on rough terrain and for evacuation by helicopter or rescue devices.

The radiotransparency ( $\approx 920$  U.H.) allows the device to be used during radiological examinations.

Decisions concerning moving and immobilising the patient, as well as the duration, method to be used and combination with other devices, shall be taken and executed by expert and trained personnel only.

The DM can be used with the following devices:

- 855.00 EVEREST BAG (using the VACUUM EVEREST KIT - 855.15)
- 870.04 LECCO 2.0 STRETCHER
- 880.01 911 CANYON STRETCHER
- 880.02 911 NET FULL STRETCHER
- 880.03 911 SHELL STRETCHER
- 890.01 ROLLY STRETCHER

### 4.2 PREPARATION

- Remove the device from the transport bag (M),
- shake it to evenly distribute the material contained in the vacuum chambers,
- check that the positioning of the straps (C and E) is suitable for the patient's height / build.

- The only numbered connections (F1) simplify the alignment of the ventral band (G),
- check that the valve (I) is closed (screwed clockwise) - fig. 3,
- place the device in the immediate vicinity of the patient.


Fig.3

## 4.3 IMMOBILISING THE PATIENT

After placing the patient on the device:

- a) connect the straps (C and E) of the same color and tension them to make the mattress take the shape of the patient (fig. 4),
- b) connect the suction pump (H) to the valve (G) (fig. 5) and suck until the necessary rigidity is obtained to ensure patient immobilization,
- c) disconnect the pump from the valve.

Avoid prolonged contact of the patient's skin with the device's materials.

Fig.4


Fig.5

Monitor the rigidity of the device, which can be:

- reduced by cuts and punctures,
- conditioned by appreciable changes in altitude (it reduces as the altitude increases and increases as the altitude decreases) and, if necessary, restore it.

## 4.4 TRANSPORTING THE PATIENT

The **VACUUM MATTRESS 4**:

- is equipped with 8 side handles for lifting and manual transport of the patient,
- if used in conjunction with the devices defined in paragraph 3.3.1, it can be raised and lowered by helicopters and life-saving devices.

## 4.5 STORAGE OF THE DEVICE

Before storing the device in the bag:

- a) unscrew the valve (fig. 6),
- b) wash as described in paragraph 5.2,
- c) fold it.


Fig.6

# CHAPTER 5

## MAINTENANCE AND REPAIR

### 5.1 GENERAL INFORMATION

The “**911 CANYON**” stretcher is made of material which is highly resistant to wear and external agents. Despite this, the use conditions require maintenance operations and even repair in particular conditions.

Maintenance and repair operations must be recorded in a special register, an example of which is indicated in chapter 10.

### 5.2 MAINTENANCE

The maintenance operations that must be carried out by the user are, in order:

- a) cleaning: after every use, wash with warm (max. 40 °C) potable water, possibly adding a neutral detergent (e.g. Marseille soap). Rinse and leave to dry in a shaded place, away from direct sources of heat,
- b) disinfection, when deemed necessary or after prolonged inactivity (more than 21 days), preceded by cleaning: wipe with a clean cloth, soaked in a water solution of sodium hypochlorite (bleach) with a 0,1% concentration (1000 ppm). When blood or other bodily fluids are present, the recommended concentration of sodium hypochlorite is 0,5% (5000ppm).

### 5.3 REPAIR

Repair shall be strictly carried out by the manufacturer.  
The user is allowed to replace the parts mentioned in paragraph 3.4.2 only, with new and genuine parts.

# CHAPTER 6

## 6 STORAGE

After cleaning, any disinfection and drying operations, keep the product and its accessories in a place that is dry (relative humidity 40-90%), cool (temperature 5-40° C) and safe (avoid UV radiation), chemically neutral (absolutely avoid salty and/or acid environments), away from sharp edges, sources of heat, humidity, corrosive substances or other possible jeopardising conditions.

**Do not store this device wet!**

# CHAPTER 7

## CHECKS INSPECTIONS AND REVISIONS

### 7.1 CHECKS

Before and after every use, check the device and ensure that:

- the pvc structure is not torn, cut, perforated or excessively worn,
- the textile parts do not show cuts, burns, residues of chemical products, excessive hair,
- markings, including labels, are readable.

### 7.2 INSPECTIONS

The device must be inspected at least once a year, starting from the date of first use, by expert personnel trained and approved by KONG S.p.A. The date of the first use and the outcome of the inspection must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.

### 7.3 REVISIONS

The device must be serviced by KONG S.p.A., or by expert personnel trained and approved by KONG S.p.A. if:

- any malfunctions are detected,
- the result of the pre and post use check or of the inspections is negative.

The outcome of the servicing must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10. The revised device is guaranteed for one year from the revision date.

# CHAPTER 8

## DURATION AND WARRANTY

### 8.1 DURATION

Carefully read the “Warnings and limitations of use” paragraph in chapter 2. The duration of this device depends on the outcome of the at least one-yearly inspection (paragraph 7.2) and possible servicing. The suspension, even if not used, has a maximum duration of 10 years from the date of manufacture indicated on the label.

**Important: any devices that have not passed the pre-use, post-use and periodic inspection must be eliminated and made unusable.**

### 8.2 DISPOSAL

Follow the rules in force in the country of use and the waste disposal procedures of the reference hospital facility for proper disposal.

### 8.3 WARRANTY

The manufacturer guarantees the compliance of the device with the regulations in force at the time of production. Warranty covering defects is limited to defects regarding raw materials and manufacturing defects. The warranty does not cover normal wear, oxidation, damage caused by non-compliant use and/or in competition, improper maintenance, transportation, preservation or storage. Warranty will be immediately deemed null and void should the device be modified or tampered with. Validity corresponds to the legal validity of the country where the device was sold, as from the date of sale, by the manufacturer. No claims shall be raised against the manufacturer once the aforementioned period has expired. Any repair or replacement request under warranty must be accompanied by a proof of purchase. Should the defect be recognised, the manufacturer undertakes to repair or, replace or refund the device at his discretion. The manufacturer shall not be held liable beyond the invoice price of the device in any case whatsoever.


### 8.4 LAW OBLIGATIONS

Professional and free time activities are often regulated by special national laws which can impose limits and/or obligations concerning the use of these devices. The user has the obligation to know and apply such laws which could provide for limits different from those indicated in this information.

# CHAPTER 9

## LABELLING AND SYMBOL

### 9.1 MEDICAL DEVICE LABEL

Any changes to the positioning of the symbols does not change the content thereof.

### 9.2 SEPARABLE COMPONENTS LABEL

Any changes to the positioning of the symbols does not change the content thereof.

### 9.3 SYMBOLS

#### Identification and traceability

	Medical Device
	Product identification code
	Unique Device Identifier
	Manufacturer's identification
	Unique serial number
 ITALY YYYY / MM	Country, year and month of production.

Data matrix UDI  
 Leggibile con le App:  
 iGepir - BarValid - Barcode Scanner

#### Information and safety

	Consult the user manual
	Compliant with Regulation (EU) 2017/745

# CHAPTER 10

## DOCUMENTS

### 10.1 INSPECTIONS AND SERVICING REGISTER

<b>VACUUM MATTRESS – Class I Medical Device</b>				
<b>REF.</b>	872.07	<b>MODEL</b>	MATTRESS VACUUM 4	
<b>UDI-DI</b>	08023577068695	<b>SN</b>		
<b>Date of first use</b>				
<b>INSPECTION AND SERVICING</b>				
<b>DATE</b>	<b>I/R</b>	<b>Description</b>	<b>Outcome</b>	<b>Supervisor</b>

# 10.2 MAINTENANCE AND REPAIR REGISTER

<b>VACUUM MATTRESS – Class I Medical Device</b>				
<b>REF.</b>	872.07	<b>MODEL</b>	MATTRESS VACUUM 4	
<b>UDI-DI</b>	08023577068695	<b>SN</b>		
<b>Date of first use</b>				
<b>MAINTENANCE AND REPAIR</b>				
<b>DATE</b>	<b>M/R</b>	<b>Description</b>	<b>Outcome</b>	<b>Supervisor</b>

# 10.3 EU DECLARATION OF CONFORMITY

The EU Declaration of Conformity is available for download at [www.kong.it](http://www.kong.it)

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[www.kong.it](http://www.kong.it)

## 11.1 APPLIED STANDARDS

- EN 597-1: 2015  
Furniture - Assessment of the ignitability of mattresses and upholstered bed bases - Part 1: Ignition source smouldering cigarette
- EN 1865-1: 2010 + A1: 2015  
Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment
- EN 62366-1: 2015 + A1: 2020  
Medical devices - Application of usability engineering to medical devices
- EN ISO 846:2019  
Plastics - Evaluation of the action of microorganisms
- EN ISO 10993-1: 2020  
Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
- EN ISO 10993-5: 2009  
Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
- EN ISO 10993-23: 2021  
Biological evaluation of medical devices - Part 23: Tests for irritation
- EN ISO 14971: 2019  
Medical devices - Application of risk management to medical devices
- EN ISO 15223-1: 2021  
Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements
- EN ISO 20417: 2021  
Medical devices - Information to be supplied by the manufacturer.

## 11.2 STANDARDS USED AS REFERENCE

- EN ISO 13485: 2016  
Medical devices - Quality management systems - Requirements for regulatory purposes
- EN 13718-1: 2014 + A1: 2020  
Medical vehicles and their equipment - Air ambulances - Part 1: Requirements for medical devices used in air ambulances.

**KONG s.p.a.**

Via XXV Aprile, 4 - (zona industriale)

I - 23804 MONTE MARENZO (LC) - ITALY

Tel +39 0341630506 - Fax +39 0341641550

**[www.kong.it](http://www.kong.it)**